

# ECO LS LINE



Available in the following versions:

	LE11L   LE11	LE16L   LE16	LE21L   LE21	LE27L   LE27
Light	with light   without light	with light   without light	with light   without light	with light   without light
Transmission	1:1	1:1	2:1	2,7:1
Label	blue	blue	green	green
Spray system	internal	internal	internal	internal
Max. speed	40.000 rpm	40.000 rpm	20.000 rpm	14.814 rpm
Connection	INTRAmatic®	INTRAmatic®	INTRAmatic®	INTRAmatic®
Suitable burs	Contra angle burs	FG burs	Contra angle burs	Contra angle burs

	LE54L   LE54	LE74L   LE74	LE148L   LE148
Light	with light   without light	with light   without light	with light   without light
Transmission	5,4:1	7,4:1	14,8:1
Label	green	green	green
Spray system	internal	internal	internal
Max. speed	7.407 rpm	5.406 rpm	2.703 rpm
Connection	INTRAmatic®	INTRAmatic®	INTRAmatic®
Suitable burs	Contra angle burs	Contra angle burs	Contra angle burs